



## ΜΙΧΑΗΛΙΔΗΣ Σ. & Ε. & ΣΙΑ Ο.Ε.

ΑΘΗΝΑ : Περικλέους 22, ΤΚ-10562, Σύνταγμα,  
▲ Τηλ. 210/3242270, Fax. 210/3255457  
e-mail : [sales@michailidis.com.gr](mailto:sales@michailidis.com.gr)

### ΣΥΣΚΕΥΗ ΑΤΜΟΥ E 503

Steam jets are a very effective method of cleaning for restoring the shine to samples, cleaning models to remove wax residues, and eliminating traces of oxidants left by welding. In the dentistry field, this method is widely used for the decontamination of metal bridges, cleaning models, washing instruments and eliminating residues left by the casting process. In general, saturated steam is also suitable for putting a gloss on paper or fabric parts, since the items always remain completely dry.

The E 5003 steam cleaner is able to produce saturated (i.e. absolutely moisture-free) steam at a pressure of 6 bar. This is essential for rapid, in-depth cleaning even of details, leaving the items for cleaning completely dry. The E 5003 is a reliable, safe, extremely compact appliance, easy to use in a thousand different applications.



#### Main features

- Safe and easy to use.
- Rapid, in -depth cleaning with saturated steam.
- Practical foot-control lever for the coming out of the steam.
- Automatic control of boiler water level.

#### Technical features

E 5003		width	depth	height
external dimensions	mm	230	350	440
max. steam pressure	bar	6		
boiler capacity	lt	4,8		
heating power	W	2000		
power supply voltage	V	220	50/60 Hz	